

Dane aktualne na dzień: 24-05-2026 21:21

Link do produktu: <https://akmed.pl/uchwyt-prosty-p-6220.html>

Uchwyt prosty



Numer katalogowy	U-200 - U-1000
------------------	-----------------------

Producent	Antar
-----------	--------------

Opis produktu

Uchwyt prosty.

Numer katalogowy:

- **U-200** uchwyt prosty **A** [mm]: 200 **B**[mm]: 95
- **U-300** uchwyt prosty **A** [mm]: 300 **B** [mm]: 95
- **U-400** uchwyt prosty **A** [mm]: 400 **B** [mm]: 95
- **U-500** uchwyt prosty **A** [mm]: 500 **B** [mm]: 95
- **U-600** uchwyt prosty **A** [mm]: 600 **B** [mm]: 95
- **U-700** uchwyt prosty **A** [mm]: 700 **B** [mm]: 95
- **U-800** uchwyt prosty **A** [mm]: 800 **B** [mm]: 95
- **U-900** uchwyt prosty **A** [mm]: 900 **B** [mm]: 95
- **U-1000** uchwyt prosty **A** [mm]: 1000 **B**[mm]: 95

Produkt posiada dodatkowe opcje:

Model: U-1000 , U-200 , U-300 , U-400 , U-500 , U-600 , U-700 , U-800 , U-900